# Dual Airbag Foot & Ankle Massager

## User Manual

Please keep this manual properly for future reference. If any abnormality occurs during use, turn off the device immediately and contact after-sales service.

### Product Introduction

This dual-airbag foot massager features a simple and stylish design and is easy to use. By rhythmically inflating and deflating the airbags, it applies cyclic compression to the ankle area, effectively soothing post-exercise muscle soreness and supporting ankle contouring. The product utilizes sequential inflation and deflation of multiple airbags to deliver even compression from the distal to the proximal part of the foot. This creates a cyclic pressure that helps dilate blood vessels, stimulate local blood circulation, invigorate tissue vitality, and improve muscle condition.

### Features

Simple operation, button controls for temperature and airbag massage.
Warm heat therapy to relax muscles.
Three massage modes to choose.
Four adjustable airbag intensity levels.

Product Parameters
Product Name: Dual Airbag Foot & Ankle Massager
Model: 101652
Battery Type: 18650
Battery Capacity: 3.7V / 3000mAh / 11.1Wh

Rated Power: 7W
Charging Current: 1A
Charging Input: 5V=1A
Charging Time: ≤3.3 hours
Usage Time: 1.5 to 5 hours
(depending on massage intensity and heating level)

### Button Description

EN

### Usage Instructions

1. Preparation: Check whether the product can be started. If not, please charge it first using the included charging cable.
2. Unfold the Product: Open the stick-on buckle and fully extend the product.
3. Wear Correctly: Wrap the product around your ankle. Ensure the main unit buttons face upward and the display screen faces downward. Adjust the tightness for comfort. Note: It can be used over clothing.
4. Power On & Adjust Intensity: Long press the power button to turn on. The default state is standby (level 0). Short press the intensity button to cycle through levels 1 to 4 (intensity increases from low to high).
5. Select Massage Mode: After powering on, short press the mode button to cycle through the three massage modes: Foot Soothing, Kneading Relaxation, and Lift Massage.
6. Activate Heating Function: After powering on, short press the heating button to activate the heating function. There are four temperature levels, indicated by LEDs from left to right (low to high): 40°C, 45°C, 50°C, 55°C (tolerance ±5°C). Heating can be used simultaneously with any massage mode.
7. Auto Shut-off: The device will automatically shut down after 15 minutes of continuous operation. Pressing any button during use will reset the timer. The device will also shut down automatically when the battery voltage drops below 2.8V to protect the battery.
8. Charging Status: The device does not operate while charging. The screen will display a scrolling number "8" during charging. When fully charged, the display shows "—" .

### Warnings

Improper operation may result in personal injury or property damage.
The following individuals should use with caution or under the guidance of a physician:
●Pregnant or postpartum women.
●Individuals with bone diseases or osteoporosis.
●Those with skin allergies, burns, or wounds.
●Cancer patients.
●Individuals with severe heart disease, hypertension, or hypotension.
●Epilepsy patients.
●Individuals with implanted electronic devices (e.g., pacemakers).

### Safety Precautions

●Start with low intensity and temperature settings. If discomfort occurs, lower the setting to a comfortable level.
●When not in use or to stop during operation, long press the power button to turn off the device.
●Do not use this product while charging; it can only be used when not charging.
●Do not use while driving.
●Do not disassemble or attempt to repair the device unless you are a qualified professional.
●If the power cord is damaged, replace it with a compatible, undamaged cord.
●Do not continue use if water or any other liquid enters the product. Do not use or store in damp environments.
●Do not use the device if it exhibits any abnormal behavior.
●Stop use immediately if you experience any discomfort.
●Store in a dry place, away from heat sources, flammable/explosive materials, and corrosive chemicals.
●Avoid contact with sharp objects to prevent punctures or damage.
●This product contains a lithium battery. Before disposal, power off the device and remove the battery. Dispose of the battery according to local regulations.
●Wireless radio frequency communication devices may interfere with the product's operation.
●Do not place or stack the product adjacent to other electronic devices.
●Use only official accessories. Non-original accessories may cause performance issues or safety hazards.

### Maintenance and Cleaning

- Main unit: Wipe the screen and casing with a soft cloth, and avoid scratches.
- Fabric parts: Do not wash with water. Use alcohol to disinfect and remove sweat stains, or wipe with a clean damp towel.

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